# CellScale biomaterials testing

## Tension Testing of Stiff Fibers

### Overview

Specimen fixation and strain measurement can be a challenge for the mechanical testing of both soft and stiff materials. This report outlines the techniques and results of testing done on a flax fibers, where the specimens are prone to breakage with clamp fixation and the strain to failure is relatively low (<2%).



### Specimen Fixation

The specimens used for this testing were single bundles of flax, each composed of a group of flax fibers. The bundle diameters were about 20µm. Previous attempts to clamp the bundles in steel jaws had resulted in specimen breakage due to mechanical damage. Attempts to clamp the bundles in soft plastic jaws had resulted in insufficient grip strength to prevent specimen pullout. The specimens were able to be successfully tested by gluing the fibers to cardboard tabs with cyanoacrylate glue. The tabs were fixed at 15mm apart prior to applying the glue in order to fix the zero specimen length.



Once mounted in the UStretch test system, a test protocol with a 400mN preload and a constant velocity stretch of 15µm/s until failure was applied. Force/displacement data and time-correlated images were collected at 15Hz during the test.

#### CELLSCALE | 11 – 564 Weber Street N., Waterloo, ON, Canada N2L 5C6 | 519-342-6870 | info@cellscale.com | www.cellscale.com

## Tension Testing of Stiff Fibers

### Results

Stress and strain were calculated for each of 3 tests using the force and displacement data collected. TO verify the strain values, the motions of the fiber were analyzed using the image analysis toolkit. The image analysis resulted in lower strains that the actuator displacementbased strain data. This is not surprising since some of the actuator displacement could have resulted in system deflection, stretching of the cardboard mounting tabs, or specimen slippage.



In all 3 cases, the specimens failed at approximately 500µm displacement and 5N of force. The resulting moduli for the 4 test (based on image analysisbased strain data were 10..8, 10.5, 7.8, and 6.3MPa.





### Conclusions

This testing shows that the UStretch is capable of measuring the stiffness of stiff fibers and that the image analysis toolkit is important for verify specimen strains.

©Copyright 2016, CellScale Biomaterials Testing

# CellScale biomaterials testing

CellScale Biomaterials Testing is the industry leader for precision biomaterial and mechanobiology test systems. Our products are being used at world-class academic and commercial organizations in over 30 countries around the globe.

Our mechanical test systems allow researchers to characterize the mechanical properties of biomaterials. Our mechanobiology technologies provide insights into the response of cells to mechanical stimulation.

CellScale's technologies are improving human health by helping researchers discover the causes of disease, improve medical treatments and devices, and advance regenerative medicine and other basic science research.

Visit our website or contact us to learn how our innovative products can help you achieve your research and development goals.

#### CELLSCALE | 11 - 564 Weber Street N., Waterloo, ON, Canada N2L 5C6 | 519-342-6870 | info@cellscale.com | www.cellscale.com